


Anti-Mre11 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P49959

Item Number: CM0024673



Introduction

Rabbit polyclonal antibody against human MRE11, validated for WB, IHC-P, and ELISA. Detects MRE11 (81 kDa) in human and mouse samples. Ideal for DNA repair and telomere research.

[Learn More](#)

Parameter	Value
Product Name	MRE11 Rabbit pAb
Remarks/Alias	ATLD; HNGS1; MRE11A; MRE11B; MRE11
Species	Human
GeneID (Human)	4361
GeneID	4361
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-302 of human MRE11 (NP_005582.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49959
Protein Weight	81kDa
Shipping	Ice bag