

Anti-Mrfap1 Polyclonal Antibody For Wb And Elisa - Q9Y605

Item Number: CM0028589

Introduction

Rabbit polyclonal antibody targeting human MRFAP1 (MORF4 family-associated protein 1), validated for Western blot and ELISA. Recombinant immunogen; molecular weight 15 kDa. Suitable for transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	MRFAP1 Rabbit pAb
Remarks/Alias	PGR1, PAM14, MRFAP1
Species	Human
GeneID (Human)	93621
GeneID	93621
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-127 of human MRFAP1 (NP_150638.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y605
Protein Weight	15kDa
Shipping	Ice bag