


Anti-Mrgprx1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Lb2

Item Number: CM0016201



Introduction

Rabbit polyclonal anti-MRGPRX1 antibody validated for WB and ELISA. Detects human MRGPRX1 (36 kDa), cross-reacts with mouse and rat. Ideal for research on itch and pain mechanisms.

[Learn More](#)

Parameter	Value
Product Name	MRGPRX1 Rabbit pAb
Remarks/Alias	GPCR, MGRG2, MRGX1, SNSR4, MRGPRX1
Species	Human
Gene ID	259249 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 263-322 of human MRGPRX1 (NP_671732.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96LB2
Protein Weight	36kDa
Shipping	Ice bag