

Anti-Mri1 Polyclonal Antibody For Wb & Elisa - Q9Bv20

Item Number: CM0014314



Introduction

Rabbit polyclonal antibody against human Mri1 (M1Pi/MTNA/MRDI) for Western blotting and ELISA. Detects a ~39 kDa band; cross-reacts with mouse and rat. Suitable for studies on methionine metabolism and cell invasion.

[Learn More](#)

Parameter	Value
Product Name	Mri1 Rabbit pAb
Remarks/Alias	M1Pi; MRDI; MTNA; Ypr118w; Mri1
Species	Human
GeneID (Human)	84245
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-369 of human Mri1 (NP_001026897.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BV20
Protein Weight	39kDa
Shipping	Ice bag