


Anti-Mrm3 Polyclonal Antibody For Wb, Elisa - Q9Hc36

Item Number: CM0027873



Introduction

Rabbit polyclonal antibody targeting human MRM3 (rRNA methyltransferase 3, mitochondrial), designed for Western blot and ELISA applications. Useful for studying mitochondrial rRNA methylation and related regulatory pathways.

[Learn More](#)

Parameter	Value
Product Name	MRM3 Rabbit pAb
Remarks/Alias	RMTL1; RNMTL1; MRM3
Species	Human
GeneID (Human)	55178
GeneID	55178
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-420 of human MRM3 (NP_060616.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HC36
Protein Weight	47kDa
Shipping	Ice bag