


Anti-Mrp1/Abcc1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P33527

Item Number: CM0004809



Introduction

Rabbit monoclonal antibody targeting human MRP1/ABCC1 (Multidrug resistance-associated protein 1), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Molecular weight: 172 kDa. Research use only.

[Learn More](#)

Parameter	Value
Product Name	MRP1/ABCC1 Rabbit mAb
Remarks/Alias	MRP; ABCC; GS-X; MRP1; ABC29; DFNA77; MRP1/ABCC1
Species	Human
GeneID (Human)	4363
GeneID	4363
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-304 of human MRP1/ABCC1 (NP_004987.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33527
Protein Weight	172kDa
Shipping	Ice bag