

Anti-Mrp1/Abcc1 Rabbit Polyclonal Antibody For Wb, Elisa - P33527

Item Number: CM0013791



Introduction

Rabbit polyclonal antibody targeting human MRP1/ABCC1, a 172 kDa multidrug resistance-associated protein. Validated for WB and ELISA. Ideal for studying drug resistance and ABC transporter function.

[Learn More](#)

Parameter	Value
Product Name	MRP1/ABCC1 Rabbit pAb
Remarks/Alias	MRP, ABCC, GS-X, MRP1, ABC29, DFNA77, MRP1/ABCC1
Species	Human
GeneID (Human)	4363
GeneID	4363
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 850-960 of human MRP1/ABCC1 (NP_004987.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P33527
Protein Weight	172kDa
Shipping	Ice bag