


Anti-Mrp1/Abcc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P33527

Item Number: CM0031092



Introduction

Rabbit polyclonal antibody targeting human MRP1/ABCC1 (172 kDa, P33527). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat MRP1. Suitable for multidrug resistance and transport studies.

[Learn More](#)

Parameter	Value
Product Name	MRP1/ABCC1 Rabbit pAb
Remarks/Alias	MRP; ABCC; GS-X; MRP1; ABC29; DFNA77; MRP1/ABCC1
Species	Human
Gene ID (Human)	4363
Gene ID	4363
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 850-960 of human MRP1/ABCC1 (NP_004987.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33527
Protein Weight	172kDa
Shipping	Ice bag