


Anti-Mrp2/Abcc2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92887

Item Number: CM0027080



Introduction

High-quality rabbit polyclonal antibody targeting human MRP2/ABCC2. Validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for studying multidrug resistance and hepatobiliary transport.

[Learn More](#)

Parameter	Value
Product Name	MRP2/ABCC2 Rabbit pAb
Remarks/Alias	DJS; MRP2; cMRP; ABC30; CMOAT; MRP2/ABCC2
Species	Human
Gene ID (Human)	1244
Gene ID	1244
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 180-315 of human MRP2/ABCC2 (NP_000383.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92887
Protein Weight	174kDa
Shipping	Ice bag