

Anti-Mrp4/Abcc4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15439

Item Number: CM0023904



Introduction

Rabbit polyclonal antibody targeting MRP4/ABCC4, a multidrug resistance transporter. Validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Ideal for studying drug resistance, cellular signaling, and transport mechanisms.

[Learn More](#)

Parameter	Specification
Product Name	MRP4/ABCC4 Rabbit pAb
Remarks/Alias	MRP4; MOATB; MOAT-B; MRP4/ABCC4
Species	Human
GeneID	10257
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1067-1325 of human MRP4/ABCC4 (NP_005836.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15439
Protein Weight	150kDa
Shipping	Ice bag