

Anti-Mrp4/Abcc4 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - O15439

Item Number: CM0005127



Introduction

KD validated rabbit monoclonal antibody targeting human MRP4/ABCC4 (multidrug resistance-associated protein 4). Ideal for WB, IF-P, IHC-P, and ELISA. Detects a 150kDa protein; immunogen covers C-terminal region (aa 1067-1325). Suitable for drug transport and resistance research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MRP4/ABCC4 Rabbit mAb
Remarks/Alias	MRP4; MOATB; MOAT-B; C4
Species	Human
GeneID (Human)	10257
GeneID	10257
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1067-1325 of human MRP4/ABCC4 (NP_005836.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O15439
Protein Weight	150kDa
Shipping	Ice bag