

Anti-Mrp4/Abcc4 Rabbit Polyclonal Antibody For Wb, Elisa - O15439

Item Number: CM0020306



Introduction

Rabbit polyclonal antibody targeting human ATP-binding cassette sub-family C member 4 (ABCC4/MRP4), validated for Western blot and ELISA. Ideal for studying drug resistance, transport, and cellular signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	MRP4/ABCC4 Rabbit pAb
Remarks/Alias	MRP4; MOATB; MOAT-B; MRP4/ABCC4
Species	Human
GeneID (Human)	10257
GeneID	10257
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1067-1325 of human MRP4/ABCC4 (NP_005836.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15439
Protein Weight	150kDa
Shipping	Ice bag