


Anti-Mrpl11 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y3B7

Item Number: CM0027484



Introduction

Rabbit polyclonal antibody targeting mitochondrial ribosomal protein MRPL11 (L11mt). Validated for WB and ELISA, reacts with human, mouse, and rat samples. Supports mitochondrial ribosome research.

[Learn More](#)

Parameter	Value
Product Name	MRPL11 Rabbit pAb
Remarks/Alias	L11MT; CGI-113; MRP-L11; MRPL11
Species	Human
GeneID	65003
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-192 of human MRPL11 (NP_057134.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3B7
Protein Weight	21kDa
Shipping	Ice bag