


Anti-Mrpl12 Rabbit Pab For Wb - P52815

Item Number: CM0025591



Introduction

Rabbit polyclonal antibody against human MRPL12, a mitochondrial ribosomal protein involved in translation and transcription activation. Validated in WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	MRPL12 Rabbit pAb
Remarks/Alias	5c5-2; L12mt; MRPL7; RPML12; MRPL7/L12; MRP-L31/34; MRPL12
Species	Human
GeneID (Human)	6182
GeneID	6182
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human MRPL12 (NP_002940.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52815
Protein Weight	21kDa
Shipping	Ice bag