


Anti-Mrpl13 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Byd1

Item Number: CM0029865



Introduction

MRPL13 rabbit polyclonal antibody validated for WB, IHC-P, ELISA with cross-reactivity to human, mouse, rat. Suitable for mitochondrial ribosome studies, protein expression analysis, and immunohistochemistry.

[Learn More](#)

Parameter	Value
Product Name	MRPL13 Rabbit pAb
Remarks/Alias	L13; L13A; L13mt; RPL13; RPML13; MRPL13
Species	Human
GeneID (Human)	28998
GeneID	28998
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human MRPL13 (NP_054797.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYD1
Protein Weight	21kDa
Shipping	Ice bag