

Anti-Mrpl13 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Byd1

Item Number: CM0019641



Introduction

High-quality anti-MRPL13 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets mitochondrial ribosomal protein L13, 20kDa. Ideal for mitochondrial function and translation research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MRPL13 Rabbit pAb
Remarks/Alias	MRPL13; L13; L13A; L13mt; RPL13; RPML13; mitochondrial ribosomal protein L13
Species	Human
GeneID (human)	28998
GeneID	28998
Immunogen	Synthetic peptide within amino acids 70-178 of human MRPL13 (NP_054797.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYD1
Protein Weight	20kDa
Shipping	Ice bag