


Anti-Mrpl15 Rabbit Polyclonal Antibody For Wb, Elisa - Q9P015

Item Number: CM0019422



Introduction

MRPL15 Rabbit Polyclonal Antibody validated for WB and ELISA, cross-reacts with human, mouse, and rat MRPL15. Ideal for mitochondrial ribosome studies and protein expression analysis.

[Learn More](#)

Parameter	Value
Product Name	MRPL15 Rabbit pAb
Remarks/Alias	MRPL15; HSPC145; L15mt; MRP-L15; MRP-L7; RPML7; mitochondrial ribosomal protein L15
Species	Human
GeneID (Human)	29088
GeneID	29088
Immunogen	Synthetic peptide, A synthetic peptide corresponding to a sequence within amino acids 197-296 of human MRPL15(NP_054894.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P015
Protein Weight	33kDa
Shipping	Ice bag