


Anti-Mrpl16 Polyclonal Antibody For Wb, Elisa - Q9Nx20

Item Number: CM0023716



Introduction

Rabbit polyclonal antibody targeting human MRPL16, a mitochondrial ribosomal protein. Validated for WB and ELISA applications in human, mouse, and rat. Ideal for mitochondrial ribosome studies.

[Learn More](#)

Parameter	Value
Product Name	MRPL16 Rabbit pAb
Remarks/Alias	L16mt; uL16m; MRP-L16; PNAS-111; MRPL16
Species	Human
GeneID (Human)	54948
GeneID	54948
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MRPL16 (NP_060310.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NX20
Protein Weight	28kDa
Shipping	Ice bag