


Anti-Mrpl18 Polyclonal Antibody For Wb, If, Elisa - Q9H0U6

Item Number: CM0019998



Introduction

MRPL18 Rabbit Polyclonal Antibody validated for Western blot (WB), IF/ICC, and ELISA.

Recognizes human and mouse MRPL18, a 21 kDa mitochondrial ribosomal protein involved in 5S rRNA import.

[Learn More](#)

Parameter	Value
Product Name	MRPL18 Rabbit pAb
Remarks/Alias	L18mt; HSPC071; MRP-L18; MRPL18
Species	Human
GeneID	29074
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human MRPL18 (NP_054880.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H0U6
Protein Weight	21kDa
Shipping	Ice bag