


Anti-Mrpl18 Rabbit Polyclonal Antibody For Wb & Elisa - Q9Cql5

Item Number: CM0024328



Introduction

Anti-Mrpl18 rabbit polyclonal antibody validated for WB and ELISA, targeting mitochondrial ribosomal protein L18 (MRP-L18). Cross-reacts with human and mouse Mrpl18, supporting mitochondrial ribosome and RNA import studies.

[Learn More](#)

Parameter	Value
Product Name	Mrpl18 Rabbit pAb
Remarks/Alias	HSPC071; MRP-L18; 1010001C05Rik; Mrpl18
Species	Human
Gene ID (Human)	29074
Gene ID	67681
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-180 of mouse Mrpl18 (NP_080586.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9CQL5
Protein Weight	21kDa
Shipping	Ice bag