

Anti-Mrpl21 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - Q7Z2W9

Item Number: CM0004294



Introduction

High-specificity rabbit monoclonal antibody targeting MRPL21, validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial protein translation and ribosome studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MRPL21 Rabbit mAb
Remarks/Alias	L21mt; bL21m; MRP-L21
Species	Human
Gene ID	219927
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 40-140 of human MRPL21 (NP_852615.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z2W9
Protein Weight	13kDa/22kDa
Shipping	Ice bag