


Anti-Mrpl22 Polyclonal Antibody For Wb, Elisa - Q9Nwu5

Item Number: CM0021848



Introduction

Rabbit polyclonal antibody against human MRPL22, suitable for Western blot and ELISA. MRPL22 is a mitochondrial ribosomal protein of the large subunit, involved in mitochondrial protein synthesis.

[Learn More](#)

Parameter	Value
Product Name	MRPL22 Rabbit pAb
Remarks/Alias	MRPL22; HSPC158; L22mt; MRP-L22; MRP-L25; RPML25; mitochondrial ribosomal protein L22
Species	Human
GeneID (Human)	29093
GeneID	29093
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-206 of human MRPL22 (NP_054899.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NWU5
Protein Weight	13kDa/14kDa/23kDa
Shipping	Ice bag