

Anti-Mrpl23 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16540

Item Number: CM0023900



Introduction

High-quality Anti-MRPL23 Rabbit Polyclonal Antibody for WB, IHC-P, IF/ICC, ELISA in Human, Mouse, Rat. Targets mitochondrial ribosomal protein L23 (~18kDa) for IVD research use.

[Learn More](#)

Parameter	Value
Product Name	MRPL23 Rabbit pAb
Remarks/Alias	RPL23; L23MRP; RPL23L; MRPL23
Species	Human
Gene ID (Human)	6150
Gene ID	6150
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-153 of human MRPL23 (NP_066957.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16540
Protein Weight	18kDa
Shipping	Ice bag