


Anti-Mrpl24 Polyclonal Antibody For Wb, Elisa - Q96A35

Item Number: CM0026160



Introduction

Rabbit polyclonal antibody targeting human mitochondrial ribosomal protein MRPL24. Validated for WB and ELISA, with cross-reactivity to mouse and rat. A tool for mitochondrial ribosome and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	MRPL24 Rabbit pAb
Remarks/Alias	L24mt; MRP-L18; MRP-L24; MRPL24
Species	Human
GeneID (human)	79590
GeneID	79590
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-216 of human MRPL24 (NP_078816.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96A35
Protein Weight	25kDa
Shipping	Ice bag