


Anti-Mrpl27 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9P0M9

Item Number: CM0021973



Introduction

KD validated anti-MRPL27 rabbit polyclonal antibody for Western blot, IHC-P, and ELISA. Detects human MRPL27, cross-reacts with mouse. Ideal for mitochondrial protein synthesis research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MRPL27 Rabbit pAb
Remarks/Alias	L27mt; MRPL27
Species	Human
GeneID	51264 (Human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 49-148 of human MRPL27 (NP_057588.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P0M9
Protein Weight	16kDa
Shipping	Ice bag