


## Anti-Mrpl28 Monoclonal Antibody For Wb, Ip, Elisa - Q13084

Item Number: CM0003680



### Introduction

High-specificity rabbit monoclonal antibody targeting MRPL28 (bL28m, p15, MAAT1). Validated for WB, IP, ELISA in human, mouse, rat. Ideal for mitochondrial ribosome and tumor antigen research.

[Learn More](#)

| Parameter        | Value                                                                                                |
|------------------|------------------------------------------------------------------------------------------------------|
| Product Name     | MRPL28 Rabbit mAb                                                                                    |
| Remarks/Alias    | p15; MAAT1; MRPL28                                                                                   |
| Species          | Human                                                                                                |
| Gene ID          | 10573                                                                                                |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MRPL28 (NP_006419.2) |
| Source           | Rabbit                                                                                               |
| Category         | Monoclonal Antibodies                                                                                |
| Application      | WB, IP, ELISA                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                    |
| SWISS            | Q13084                                                                                               |
| Protein Weight   | 30kDa                                                                                                |
| Shipping         | Ice bag                                                                                              |