


Anti-Mrpl28 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q13084

Item Number: CM0028710



Introduction

High-quality rabbit anti-MRPL28 polyclonal antibody validated for WB, IHC-P, IF/ICC, IP, and ELISA applications in human and rat samples. Ideal for mitochondrial ribosome and melanoma-associated antigen studies.

[Learn More](#)

Parameter	Specification
Product Name	MRPL28 Rabbit pAb
Remarks/Alias	p15; MAAT1; MRPL28
Species	Human
GeneID (Human)	10573
GeneID	10573
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-256 of human MRPL28 (NP_006419.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13084
Protein Weight	30kDa
Shipping	Ice bag