


Anti-Mrpl30 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Tcc3

Item Number: CM0024307



Introduction

Affinity-purified rabbit polyclonal antibody against human MRPL30, validated for Western blot, IHC-P, and ELISA. Cross-reacts with rat. Suitable for mitochondrial ribosome and translation studies.

[Learn More](#)

Parameter	Value
Product Name	MRPL30 Rabbit pAb
Remarks/Alias	L28MT; L30MT; MRPL28; RPML28; MRP-L28; MRP-L30; MRPL28M; MRPL30
Species	Human
Gene ID	51263
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-161 of human MRPL30 (NP_660213.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8TCC3
Protein Weight	19kDa
Shipping	Ice bag