

Anti-Mrpl32 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Byc8

Item Number: CM0016931



Introduction

Rabbit polyclonal antibody targeting human MRPL32, a mitochondrial ribosomal protein essential for mitoribosomal protein synthesis. Validated for WB, IF/ICC, ELISA; cross-reacts with mouse. Suitable for mitochondrial research and diagnostics.

[Learn More](#)

Parameter	Value
Product Name	MRPL32 Rabbit pAb
Remarks/Alias	L32mt; HSPC283; MRP-L32; bMRP-59b; MRPL32
Species	Human
Gene ID (Human)	64983
Gene ID	64983
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 19-188 of human MRPL32 (NP_114109.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BYC8
Protein Weight	21kDa
Shipping	Ice bag