


Anti-Mrpl33 Rabbit Polyclonal Antibody For Ihc-P, Elisa - 075394

Item Number: CM0020174



Introduction

MRPL33 Rabbit Polyclonal Antibody validated for IHC-P and ELISA. Targets human mitochondrial ribosomal protein L33 (bL33m), a component of the 39S large subunit. Suitable for mitochondrial translation and disease research.

[Learn More](#)

Parameter	Specification
Product Name	MRPL33 Rabbit pAb
Remarks/Alias	L33mt; bL33m; C2orf1; RPL33L; MRP-L33; MRPL33
Species	Human
Gene ID (Human)	9553
Gene ID	9553
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-65 of human MRPL33 (NP_004882.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	O75394
Protein Weight	7kDa
Shipping	Ice bag