

Anti-Mrpl38 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Dv4

Item Number: CM0021478



Introduction

High-quality anti-MRPL38 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Detects human, mouse, and rat MRPL38, a mitochondrial large ribosomal subunit protein. Ideal for mitochondrial translation and protein expression studies.

[Learn More](#)

Parameter	Specification
Product Name	MRPL38 Rabbit pAb
Remarks/Alias	L38MT; RPML3; MRP-L3; HSPC262; MRP-L38; MRPL38
Species	Human
Gene ID (Human)	64978
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96DV4
Protein Weight	45kDa
Shipping	Ice bag