


Anti-Mrpl39 Polyclonal Antibody For Wb - Q9Nyk5

Item Number: CM0014379



Introduction

Rabbit polyclonal antibody specific to MRPL39, a mitochondrial large ribosomal subunit protein. Validated for WB and ELISA analysis in human, mouse, and rat samples. Ideal for mitochondrial translation studies.

[Learn More](#)

Parameter	Value
Product Name	MRPL39 Rabbit pAb
Remarks/Alias	L5mt; L39mt; MRPL5; RPML5; MRP-L5; PRED22; PRED66; MSTP003; C21orf92; MRPL39
Species	Human
Gene ID (Human)	54148
Gene ID	54148
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-338 of human MRPL39 (NP_059142.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYK5
Protein Weight	39kDa
Shipping	Ice bag