


Anti-Mrpl40 Polyclonal Antibody For Wb, Elisa - Q9Nq50

Item Number: CM0021409



Introduction

High-quality rabbit polyclonal antibody against human MRPL40 (mL40, mitochondrial ribosomal protein L40). Suitable for WB, ELISA. Detects endogenous MRPL40 at ~24 kDa. Validated in human samples for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	MRPL40 Rabbit pAb
Remarks/Alias	URIM; L40mt; NLVCF; MRPL22; MRP-L22; MRP-L40; MRPL40
Species	Human
GeneID (Human)	64976
GeneID	64976
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-206 of human MRPL40 (NP_003767.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NQ50
Protein Weight	24kDa
Shipping	Ice bag