


Anti-Mrpl42 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y6G3

Item Number: CM0018455



Introduction

Anti-MRPL42 Rabbit Polyclonal Antibody validated for Western Blot and ELISA. Reacts with Human, Mouse. Recognizes mitochondrial ribosomal protein L42 (17kDa). Ideal for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	MRPL42 Rabbit pAb
Remarks/Alias	L31MT; L42MT; S32MT; MRPL31; MRPS32; PTD007; RPML31; HSPC204; MRP-L31; MRP-L42; MRP-S32; MRPL42
Species	Human
GeneID (Human)	28977
GeneID	28977
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human MRPL42 (NP_054769.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y6G3
Protein Weight	17kDa
Shipping	Ice bag