


Anti-Mrpl43 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8N983

Item Number: CM0019073



Introduction

Rabbit polyclonal antibody targeting human MRPL43, a mitochondrial ribosomal protein. Suitable for Western blot, immunofluorescence, and ELISA; cross-reacts with human, mouse, and rat. Ideal for mitochondrial translation studies.

[Learn More](#)

Parameter	Value
Product Name	MRPL43 Rabbit pAb
Remarks/Alias	L43mt; MRP-L43; bMRP36a; MRPL43
Species	Human
GeneID (human)	84545
GeneID	84545
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-155 of human MRPL43 (NP_789762.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N983
Protein Weight	23kDa
Shipping	Ice bag