


Anti-Mrpl46 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Eqi8

Item Number: CM0023475



Introduction

Mrpl46 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA in human, mouse, and rat. Targets mitochondrial large ribosomal subunit protein mL46 (39S ribosomal protein L46), a component of the mitochondrial ribosome essential for mitochondrial translation.

[Learn More](#)

Parameter	Value
Product Name	Mrpl46 Rabbit pAb
Remarks/Alias	LIECG2; P2ECSL; C15orf4; 3110052F15Rik; Mrpl46
Species	Human
Gene ID (Human)	26589
Gene ID	67308
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of mouse Mrpl46 (NP_075820.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9EQI8
Protein Weight	32kDa
Shipping	Ice bag