

Anti-Mrpl46 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H2W6

Item Number: CM0018557



Introduction

Rabbit polyclonal antibody targeting human MRPL46, a mitochondrial ribosomal protein mL46. Validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for mitochondrial translation studies.

[Learn More](#)

Parameter	Value
Product Name	MRPL46 Rabbit pAb
Remarks/Alias	LIECG2; P2ECSL; C15orf4; MRPL46
Species	Human
Gene ID	26589
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-279 of human MRPL46 (NP_071446.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H2W6
Protein Weight	32kDa
Shipping	Ice bag