


Anti-Mrpl50 Polyclonal Antibody For Wb, Elisa - Q8N5N7

Item Number: CM0024394



Introduction

Rabbit polyclonal antibody against human MRPL50 (mL50, MRP-L50), a mitochondrial ribosomal protein. Validated for WB and ELISA, cross-reacts with human. Useful for mitochondrial translation studies.

[Learn More](#)

Parameter	Value
Product Name	MRPL50 Rabbit pAb
Remarks/Alias	mL50; MRP-L50
Species	Human
GeneID (human)	54534
GeneID	54534
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MRPL50(NP_061924.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8N5N7
Protein Weight	18kDa
Shipping	Ice bag