


## Anti-Mrpl51 Polyclonal Antibody For Wb, Elisa - Q4U2R6

Item Number: CM0016629



### Introduction

Rabbit polyclonal antibody against human mitochondrial ribosomal protein MRPL51, validated for Western blotting and ELISA. Cross-reacts with mouse and rat MRPL51. Ideal for mitochondrial ribosome studies.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | MRPL51 Rabbit pAb                                                                                                 |
| Remarks/Alias    | CDA09, MRP64, bMRP64, HSPC241, MRPL51                                                                             |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 51258                                                                                                             |
| Gene ID          | 51258                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-128 of human MRPL51 (NP_057581.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Mouse, Rat                                                                                                        |
| SWISS            | Q4U2R6                                                                                                            |
| Protein Weight   | 15kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |