


## Anti-Mrpl54 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6P161

Item Number: CM0014748



### Introduction

MRPL54 Rabbit pAb specifically detects human mitochondrial large ribosomal subunit protein mL54 (L54mt). Validated for WB, IHC-P, and ELISA applications. Suitable for mitochondrial translation and ribosome research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | MRPL54 Rabbit pAb                                                                                                 |
| Remarks/Alias    | L54mt; MRP-L54; MRPL54                                                                                            |
| Species          | Human                                                                                                             |
| GeneID (human)   | 116541                                                                                                            |
| GeneID           | 116541                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-138 of human MRPL54 (NP_758455.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IHC-P, ELISA                                                                                                  |
| Cross Reactivity | Human                                                                                                             |
| SWISS            | Q6P161                                                                                                            |
| Protein Weight   | 16kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |