


## Anti-Mrpl55 Polyclonal Antibody For Wb And Elisa - Q9Cz83

Item Number: CM0024137



### Introduction

High-quality rabbit polyclonal antibody against mitochondrial ribosomal protein Mrpl55. Suitable for WB and ELISA with cross-reactivity to human and mouse. Ideal for mitochondrial research and protein synthesis studies.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Mrpl55 Rabbit pAb                                                                                                     |
| Remarks/Alias    | L55mt; MRP-L55; 2810038N09Rik; Mrpl55                                                                                 |
| Species          | Human                                                                                                                 |
| Gene ID (Human)  | 128308                                                                                                                |
| Gene ID          | 67212                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 33-127 of mouse Mrpl55 (NP_001289265.1) |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse                                                                                                          |
| SWISS            | Q9CZ83                                                                                                                |
| Protein Weight   | 15kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |