


Anti-Mrpl55 Polyclonal Antibody For Wb, Elisa - Q7Z7F7

Item Number: CM0019525



Introduction

High-quality rabbit polyclonal antibody against human MRPL55, a mitochondrial ribosomal protein. Validated for WB and ELISA, ideal for studying mitochondrial translation and protein synthesis.

[Learn More](#)

Parameter	Value
Product Name	MRPL55 Rabbit pAb
Remarks/Alias	L55nt; MRP-L55; AAVG5835; PRO19675; MRPL55
Species	Human
Gene ID (Human)	128308
Gene ID	128308
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-128 of human MRPL55 (NP_001308213.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q7Z7F7
Protein Weight	15kDa
Shipping	Ice bag