


Anti-Mrpl9 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Byd2

Item Number: CM0018832



Introduction

KD Validated rabbit polyclonal antibody against mitochondrial ribosomal protein MRPL9, suitable for WB, IF/ICC, and ELISA. Detects endogenous MRPL9 in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MRPL9 Rabbit pAb
Remarks/Alias	L9mt; MRPL9
Species	Human
GeneID (Human)	65005
GeneID	65005
Immunogen	Recombinant protein corresponding to a sequence within amino acids 53-267 of human MRPL9 (NP_113608.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYD2
Protein Weight	30kDa
Shipping	Ice bag