


Anti-Mrps11 Monoclonal Antibody For Wb, If-P, Elisa - P82912

Item Number: CM0005185



Introduction

Rabbit monoclonal antibody targeting human MRPS11, a mitochondrial ribosomal protein. Validated for WB, IF-P, ELISA, and cross-reactive with human, mouse, and rat. Ideal for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	MRPS11 Rabbit mAb
Remarks/Alias	HCC-2; S11mt; MRP-S11
Species	Human
GeneID (human)	64963
GeneID	64963
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 95-194 of human MRPS11(NP_001308899.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P82912
Protein Weight	21kDa
Shipping	Ice bag