


Anti-Mrps11 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P82912

Item Number: CM0023147



Introduction

Rabbit polyclonal antibody against MRPS11, validated for Western blot, immunofluorescence, and ELISA. Detects mitochondrial ribosomal protein S11 in human, mouse, and rat. Useful for mitochondrial function and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MRPS11 Rabbit pAb
Remarks/Alias	MRPS11; HCC-2; MRP-S11; S11mt; mitochondrial ribosomal protein S11
Species	Human
GeneID (Human)	64963
GeneID	64963
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 95-194 of human MRPS11 (NP_001308899.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P82912
Protein Weight	16kDa/20kDa
Shipping	Ice bag