


Anti-Mrps12 Rabbit Polyclonal Antibody For Wb, Elisa - O15235

Item Number: CM0017880



Introduction

MRPS12 rabbit polyclonal antibody validated for Western blot and ELISA. Targets human mitochondrial ribosomal protein S12 (O15235), an essential small-subunit component for mitochondrial protein synthesis. Ideal for mitochondrial dysfunction studies.

[Learn More](#)

Parameter	Value
Product Name	MRPS12 Rabbit pAb
Remarks/Alias	RPS12; RPMS12; RPSM12; MPR-S12; MT-RPS12; MRPS12
Species	Human
GeneID (human)	6183
GeneID	6183
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-138 of human MRPS12 (NP_203527.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15235
Protein Weight	15kDa
Shipping	Ice bag