


Anti-Mrps14 Rabbit Polyclonal Antibody For Wb, Elisa - O60783

Item Number: CM0015358



Introduction

Rabbit polyclonal antibody targeting human MRPS14 (uS14m), a mitochondrial ribosomal protein. Suitable for WB and ELISA, with cross-reactivity to mouse and rat. Useful for mitochondrial translation and disease studies.

[Learn More](#)

Parameter	Value
Product Name	MRPS14 Rabbit pAb
Remarks/Alias	S14mt; COXPD38; MRP-S14; HSMRPS14; DJ262D12.2; MRPS14
Species	Human
GeneID	63931
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-128 of human MRPS14 (NP_071383.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O60783
Protein Weight	15kDa
Shipping	Ice bag