


## Anti-Mrps15 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P82914

Item Number: CM0006187



### Introduction

Anti-MRPS15 rabbit monoclonal antibody for WB, IHC-P, IF/ICC, ELISA. Recognizes human mitochondrial ribosomal protein S15 (MRP-S15, uS15m), a component of the 28S small ribosomal subunit. Suitable for diagnostic and research applications.

[Learn More](#)

| Parameter        | Value                                                                                                |
|------------------|------------------------------------------------------------------------------------------------------|
| Product Name     | MRPS15 Rabbit mAb                                                                                    |
| Remarks/Alias    | DC37, S15mt, RPMS15, MPR-S15, MRPS15                                                                 |
| Species          | Human                                                                                                |
| GeneID           | 64960                                                                                                |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 158-257 of human MRPS15 (P82914). |
| Source           | Rabbit                                                                                               |
| Category         | Monoclonal Antibodies                                                                                |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                             |
| Cross Reactivity | Human                                                                                                |
| SWISS            | P82914                                                                                               |
| Protein Weight   | 30kDa                                                                                                |
| Shipping         | Ice bag                                                                                              |