


Anti-Mrps18A Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nvs2

Item Number: CM0022655



Introduction

High-quality anti-MRPS18A (mL66) rabbit polyclonal antibody for Western blot and ELISA. Reacts with human and mouse. Ideal for mitochondrial ribosomal protein function research. 22 kDa. Q9NVS2.

[Learn More](#)

Parameter	Value
Product Name	MRPS18A Rabbit pAb
Remarks/Alias	mL66; S18bmt; S18mt-a; MRPS18-3; HumanS18b; MRP-S18-3; MRP-S18-a; MRPS18A
Species	Human
GeneID (Human)	55168
GeneID	55168
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 97-196 of human MRPS18A (NP_060605.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NVS2
Protein Weight	22kDa
Shipping	Ice bag