

Anti-Mrps18B Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y676

Item Number: CM0018119



Introduction

High-quality rabbit polyclonal antibody against human MRPS18B (mitochondrial ribosomal protein mS40). Suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Useful for mitochondrial translation and ribosome research.

[Learn More](#)

Parameter	Details
Product Name	MRPS18B Rabbit pAb
Remarks/Alias	PTD017; S18amt; C6orf14; HSPC183; MRPS18-2; HumanS18a; MRP-S18-2; MRPS18B
Species	Human
GeneID (Human)	28973
GeneID	28973
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-258 of human MRPS18B (NP_054765.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y676
Protein Weight	29kDa
Shipping	Ice bag